

CCP3A-2217-EU Care Call Pager Unit

Lightweight & discreet, this fully portable vibrating pager unit with internal rechargeable battery is worn by the user on the belt, in a pocket or around the neck using the cord supplied.

The pager can be switched between audio & vibration alert.

The units internal battery is recharged using the Care Call trickle charger and will provide 24 hour response to notification from the full range of Care Call monitors.



Alert Options include

Person to Person Key Fob	Magnetic Door Monitor
Chair leaving alarm	PIR Movement Monitor
Fall Alarm	Bed Leaving Alarm
Audio Monitor	Low Monitor battery

The unique Surekey™ feature, together with a strong vibration and night time vibrating pillow pad, makes the pager the ideal solution to aid the Carer by giving a quick and positive notification of the activation of any of the monitors coded to that pager.

Where additional monitors of the same type are required, for example when monitoring several doors or individuals in different rooms, it is simple to designate different keys on the Pager for each monitor. Up to 12 individual monitors can be programmed to illuminate different keys, giving the Carer plenty of flexibility when setting up the system.



- Lightweight and discreet
- Clear backlit display
- Vibration or Audio alert
- Integral battery level display
- Last event recall
- Long range option extends the open air range from 100m to around 1000m
- Monitor low battery notification
- 1 year warranty

New Audio Alert

The pager will now beep or discreetly vibrate when a signal is received.




Audio

or



Vibration

Operating Range



CCS3A-2217-EU SignWave Portable Sounder Flasher

Simple to operate and ergonomically designed for use by Carers with varying requirements, the portable SignWave Receiver uses an integral strobe, adjustable sounder, and the numbered display to notify the user to an event.

For use both in the home and Care facilities, the SignWave is designed to give prompt notification of any of the monitored events. With the optional vibrating pillow pad connected to the SignWave charger, the SignWave gives 24hour warning of monitored events, even for the heaviest of sleepers.



Alert Options include

Person to Person Key Fob	Magnetic Door Monitor	Chair leaving alarm
PIR Movement Monitor	Fall Alarm	Bed Leaving Alarm
Audio Monitor	Low Monitor battery	

Where a number of monitors are required, for example when monitoring several doors or individuals in different rooms, it is simple to designate different numbers on the SignWave for each monitor.

Up to 12 individual monitors can be programmed to illuminate different keys on the SignWave, giving the Carer plenty of flexibility when setting up the system.

+6

6

5

4

3

2

1

- Simple event display
- Ergonomic design making it easy to pick up and dock
- Large cancel button
- Volume control
- Options to cancel strobe or Sound
- Pillow Pad option
- Integral power cord and Pillow Pad locking mechanism
- Mains failure backup
- Long range option extends the open air range from 100m to approximately 1000m
- 1 year warranty

Operating Range




CCF3A-2219-EU Care Call Person to Person Key Fob

CCF3A-2219-EU-FALL Care Call Person to Person Key Fob (with shock detector)

The Mini Person to Person Key Fob can be used to page a carer who is in another part of the home or care facility: ideal when caring for someone who is either ill, has mobility problems or is prone to falling.


Powered by a 6V alkaline battery, the Key Fob when operated will send a signal to a Pager or SignWave Portable Receiver. The Key Fob will also indicate to the Receiver automatically when the battery needs replacing.



- Simple to use
- Neck cord incorporating a safety break
- Automatic low battery signal to Pager or SignWave
- Shock detector option to alert a carer to a fall
- 1 year warranty
- Operating range in open air approximately 1000m



Operating Range



CCM4A-2214-EU Bed Leaving Monitor

CCS-PAD Under-Carpet Pressure Pad

This compact battery operated wireless Monitor is designed to be used with an under carpet pressure pad. When a person steps onto the Pressure Pad, the monitor sends a signal to the Pager or Portable SignWave.

The Mini Monitor is ideal for use when caring for someone who may wander, or need assistance when getting out of bed.


The Monitor is also ideal in any situation where a person may need to be alerted to someone entering or leaving a room either in the home or a care facility.



- Simple to install
- Discreet
- Automatic low battery signal to Pager or SignWave
- 1 year warranty
- Simple screw fixings included
- Auto battery voltage indication on start-up
- Operating range in open air approximately 1000m



Operating Range



CCC4A-2212-EU Chair Leaving Mini Monitor

CCS-CHAIR-PAD Chair Leaving Pad

This compact, battery operated wireless Chair Leaving Monitor is both simple to fit and discreet in operation. Powered by a 9V alkaline battery, the Monitor sends a signal to a Pager or SignWave Portable Receiver when a person releases pressure from the Chair Leaving pad.

The Chair leaving Monitor is ideal for use when caring for someone who may wander or have difficulty when trying to stand.


When a signal is received from the Monitor it will be indicate a designated key on the Receiver. If a second Monitor is required, another key can be designated.



- Simple to install
- Discreet
- Automatic low battery signal to Pager or SignWave
- 1 year warranty
- Simple screw fixings included
- Auto battery voltage indication on start-up
- Operating range in open air approximately 1000m



Operating Range



CCD4A-2220-EU Care Call Magnetic Door Monitor


This compact battery operated wireless door monitor is simple to fit and discreet. Powered by a 9V battery, the unit will send a signal to the Pager or Portable SignWave each time the magnet is moved away from the Monitor.

The Magnetic Door Monitor is ideal for use when caring for someone who may wander, or in any situation where a person may need to be alerted to the opening of a door or window.


The Monitor will also alert the user when the internal battery needs replacing.



- Simple to install
- Discreet
- Automatic low battery signal to Pager or SignWave
- 1 year warranty
- Simple screw fixings included
- Auto battery voltage indication on start-up
- Operating range in open air approximately 1000m



Operating Range



CCM3A-2212-EU Audio Monitor CCM3A-2214-EU-LEAD Audio Monitor with Lead

The Audio Monitor, which can be battery or mains powered, allows the user to monitor a room or area by responding to pre-determined levels of sound. When a sound of the correct volume and duration is picked up by the monitor, it sends a signal to the Pager or SignWave portable Receiver.

By using the trigger delay keys and sensitivity control to find the optimum setting, the Carer can easily determine the level and frequency of sound a cared-for person may make, in the same way a parent would monitor a child with a baby alarm. Visual confirmation of the Audio Monitoring function can be selected if required and is ideal for user setup.

The optional Monitor with lead allows for easy connection to additional devices such as pressure pads and switching devices.



Silent Alert International DO NOT recommended the use of audio monitoring of smoke Alarms.



- Battery or Mains Adaptor Options
- Sensitivity Control
- Three trigger delays: - Short, Medium and Long
- Automatic battery voltage indication on start up
- Automatic low battery signal to Pager or SignWave
- Audio monitoring display
- Can be set to light any symbol on Pager or SignWave
- 1 year warranty
- Operating range in open air approximately 1000m

Operating Range



CCM3A-2212-EU-PIR Movement Monitor with audio sensing

This compact mains powered wireless Movement Monitor is simple to fit and easy to install. Using an infrared sensor, the Monitor sends a signal to the Pager or Portable SignWave receiver whenever a person passes through the sensors field of view.

Ideal for use when caring for someone who may wander or young child, the Monitor can be angled to cover a specific room or staircase, but allow's free passage to other rooms without activating the Pager or SignWave Receiver.

The Monitor can also be located low down at a bedside to act as a bed leaving alarm where a pressure pad would be inappropriate.


The unit features 3 non return delays of 5,10 and 15 minutes if the unit is triggered. For example if a person goes to the toilet and returns within the designated time, the unit will not send a signal to the carer.

The Movement Monitor also has an Audio sensor with variable sensitivity and time delays for use when monitoring sounds within a room such as a child crying or an elderly person calling for assistance.



- Simple to install
- Adjustable directional sensor
- Audio Sensor with volume sensitivity and time delay
- Operating range in open air approximately 1000m
- 1 year warranty
- Simple screw fixings included
- Optional low voltage extension lead available

Operating Range



JELLY-BEAN Jelly Bean Twist Switch

CCM4A-2214-EU Mini Monitor with Socket

The Jelly Bean switch enables people with poor dexterity or restricted movement to be able to raise an alert to their Carer by simply using a Jelly Bean switch connected to a standard Care Call Mini Monitor with socket. With a diameter of 64mm, the Jelly Bean offers a much larger contact area than the Person to Person Key Fob or Mini Monitor buttons, but is not too large so as to be out of place next to a bedside or chair.

There are four different coloured tops included with each Jelly Bean switch, which allows for additional switches to be linked into separate Mini Monitors to provide the option of multiple alerts.

A Carer could, for instance, use a red top on one switch to call for assistance with toileting needs and perhaps a second switch with a blue top to request a drink or for help with getting out of bed. Each of these alerts can be set to a different indicator on the Carer's Pager or SignWave Portable Receiver thereby helping the Carer to save valuable time when going to assist the person in care.




- Supplied with Red, Yellow, Blue and Green switch tops as standard.
- Extra colours available.
- Designed for general input switch applications
- Injection moulded to produce a very robust product
- Easy to mount with standard flange base
- Auditory & Tactile feedback



Other switch options are available please contact our sales department for further information.

Operating Range



CCM4A-2214-EU Mini Monitor -Special CCM3A-2214-EU Universal Monitor-Special

Silent Alert Section 7 Special Monitors have been adapted for use with specific third party products to enable Carers to be alerted to the activation of monitors or alarms of non Silent Alert International equipment and of events that they would otherwise not be able to monitor.

Developed by our Electronics' Team, Section 7 Monitors are designed to interface with an extensive range of equipment by the addition of bespoke connectors and software programmes that will enable an alert to be sent to a Pager or SignWave Portable Receiver.

The Monitors can be set-up to inform the Carer of events such as the activation of an Epilepsy monitor, Bed occupancy monitor and Neonatal Respiration monitor alarms. Many specialised alerting applications can be accommodated.



- Epilepsy alarms
- Bed Occupancy Alarms
- Respiratory equipment alarms
- Enuresis Alarms
- Automatic low battery signal to Pager or SignWave
- 1 year warranty
- Operating range in open air approximately 1000m

