

CCC4A-2212-EU Chair Leaving Mini Monitor

CCS-CHAIR-PAD Chair Leaving Pad

This compact, battery operated wireless Chair Leaving Monitor is both simple to fit and discreet in operation. Powered by a 9V alkaline battery, the Monitor sends a signal to a Pager or SignWave Portable Receiver when a person releases pressure from the Chair Leaving pad.

The Chair leaving Monitor is ideal for use when caring for someone who may wander or have difficulty when trying to stand.

When a signal is received from the Monitor it will be indicate a designated key on the Receiver. If a second Monitor is required, another key can be designated.



- Simple to install
- Discreet
- Automatic low battery signal to Pager or SignWave
- 1 year warranty
- Simple screw fixings included
- Auto battery voltage indication on start-up
- Operating range in open air approximately 1000m



Operating Range

